

cedis Hygiene and
Care Products for
everydayHEARING



egger



While wearing hearing systems, earmolds, hearing protection or similar products, the build-up of surface contamination or the accumulation of moisture is not unusual. However, it would be unusual not to act effectively against them.

For hygienic cleanliness and the preservation of function and value the right auxiliary product is necessary. Just as natural as hand washing is for personal hygiene.

egger offers the right, proven cedis brand product for almost all wishes and needs.

For more than three decades, today and tomorrow.




The fresh-colored cedis logo demands “Yes to Hygiene and Care for everydayHEARING!” Ideally, with the cleaning, function- and value-retaining cedis brand products by egger.

The right product for your needs.




In-the-ear hearing systems (ITE)
Electronic in-ear hearing protection



Insert earphones




Behind-the-ear hearing systems (BTE)




Custom-made earmolds
Custom-made hearing protection
Custom-made ear protection
Reusable earplugs

The practical starter sets, e. g.

cedis cleaning and drying set BTE




cedis cleaning and drying set ITE



Please contact your hearing care
professional for information on other
cedis products and sets.

cedis e300.DUS electrical hygiene station

- electrical, gentle drying depending on model, with miniature fan
- hygienic cleaning with UV-C light
- safe storage for hearing systems and earmolds
- automatic switch-off at program end




Please take additional information
from our detailed product leaflet.
Please contact your hearing care
professional.

cedis e150.DUS

electrical “good + affordable”
Hygiene Station

- hygienic cleaning with two opposing UV-C light sources
- electrical, gentle drying
- moderate heat generation, especially for the drying of rechargeable hearing systems




Please take additional information
from our detailed product leaflet.
Please contact your hearing care
professional.

cedis e100.DS

electrical drying station


- sufficient room for drying and storing, e. g. hearing systems with earmolds, Cochlea Implant BTE processors, earphones etc.
- the air is evenly warmed via an electrical heating control and a miniature fan
- automatic switch-off after the end of the drying process



Please take additional information
from our detailed product leaflet.
Please contact your hearing care
professional.


cedis disinfectant spray

- for simple, quick and thorough cleaning
- with either atomizer, brush or as pocket model
- for applying onto a disposable wipe and then cleaning the surface



cedis disinfectant wipes

- ready-to-use, moist wipes for surface cleaning
- practical – available in various sizes
- single wipes, pocket dispenser, large dispenser




cedis cleansing tablets

- dissolve in water
- remove even persistent deposits and contamination on the surface
- for use with cedis cleaning container and air blower for drying (available separately and as set)




cedis drying capsules

- remove accumulated perspiration or condensation
- gentle to the sensible electronics of the hearing systems
- drying process is best done over night
- for use with an airtight cedis drying container or drying bag (available separately and as set)



cedis care gel

- skin care gel with feel-good effect
- with Aloe Vera to prevent irritations
- facilitates insertion



A L O E V E R A



cedis®



Worth more than a 1000 words:
cedis 3D animation films
www.egger-info.com

796451105-2019-12

We are pleased to advise you:

is a registered trademark of
egger Otoplastik + Labortechnik GmbH



egger Otoplastik + Labortechnik GmbH · Aybühlweg 59 · 87439 Kempten/Germany
www.egger.online